Jumper Medical

JPD-ES100

Manuál

Obsah balení

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Manual

Contents of the package

Electrostimulator, two electrode sizes, charging cable, manual

Recharge

The device recharges using the supplied microUSB cable, which you connect to the charger at the other end.

Fixing electrodes

The electrodes are fixed using a magnet on the lower part of both parts of the electrostimulator. Choose the size of the electrodes according to the size of the area of the body to which you want to apply the treatment.

Charm

You turn on the electrostimulator using a button located on one of the parts of the electrostimulator.

This is the device's only charm. In addition to switching on, it is also used for switching off and switching individual treatment methods.

Connect to a mobile app

Turn on the electrostimulator and start jumperhealth. Here, start the measurement using an electrostimulator and the device will automatically connect to your mobile phone. The link works on the basis of Bluetooth. In the History tab, you can see the history of treatment.

Use

Remove the foil from the electrodes and place them in the desired place on the body. Turn on the device and pair it with the app. Note the device works without linking to the application, but you cannot control the duration of treatment, which is set by default to 15 minutes.

Set the duration of treatment in the app. At the same time, you can see what methods of treatment you have set up (it switches with the button on the device).

The intensity of treatment is always set to the mildest degree at the beginning. We strongly recommend adding intensity carefully and only after making sure that the intensity is weak that you add another degree.

Note that if you change its method during treatment, the intensity automatically switches to the lowest value.

After treatment, the device will automatically switch off.

Diode on the device

Green color – the device is charged and in operation

Flashing green – the device will turn off

Flashing Orange in no time – low battery capacity / electrode does not hold on the body

Orange color – charging

Safety instructions

Do not submerge the device in water. Clean the device with a gently damp cloth. Do not use the device if it shows signs of any damage. Do not leave the device in direct sunlight. After using the electrodes, place the protective film back on them.

Disposal

The product contains batteries and/or recyclable electrical waste. In order to protect the environment, do not dispose of the product with ordinary waste, but hand it over for recycling to a collection point of electrical waste.

Declaration of conformity

The importer, Satomar s.r.o (Karlova 37, 614 00 Brno), declares that the product meets all the essential requirements of eu directives that apply to it.

 The full text of the Declaration of Conformity can be downloaded <http://www.satomar.cz>.

Copyright © 2021 Satomar, s. r. o. All rights reserved.
Appearance and specifications may be changed without notice. Print errors reserved.

The most up-to-date version of the manual can always be found www.mobilnipodpora.cz.

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